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### Simultaneous Determination of Oxygen and Carbon Dioxide Abundance Ratio in Package by Headspace Gas Analysis Technique

CHEN Xin<sup>1</sup>, WU Xiong-jie<sup>2</sup>, CHEN Xiao-nan<sup>3</sup>, ZHOU Wei-fang<sup>1</sup>

(1. Labthink Instruments, Co., Ltd., Jinan 250031, China;

2. Anhui Supervision and Inspection Center for Package Printing Product, Tongcheng, 231400, China;

3. Jinan Food Quality Supervision and Inspection Center, Jinan 250002, China)

**Abstract:** Oxygen and carbon dioxide abundance ratio of sample with different packaged form were tested with selective test of sampling time and analysis time. Stability and application scope of headspace gas analysis technique for simultaneous determination of oxygen and carbon dioxide abundance ratio in food or pharmaceutical package was verified. The purpose was to provide reference for food, drug manufacturing enterprise or research

institute on adjusting appropriate proportion of gases in package.

**Key words:** headspace gas analysis; simultaneous determination; oxygen and carbon dioxide abundance ratio; sampling time; analysis time of sample gas

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Tab.1 Choice of sampling or analysis time in experiment process

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Tab. 2 Comparison of headspace gas data of sample with different packaged forms

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Tab. 3 Comparison of experimental data repeatability of sample with different packaged forms

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WANG Xing-dong,ZHAO Jiang.Packing Efficiency and Correlation detection of Nitrogen[J].Packaging Engineering, 2007,28(9):213-214.

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E-mail: [marketing@labthink.cn](mailto:marketing@labthink.cn)

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FU Xiao-hua.Study on Mixing and Detecting of Gas Proportion in MAP System[D].Harbin: Harbin University of Commerce,2011.

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ZHAO Su-fen.Application of Modified Atmosphere Packing on Chilled Fresh Pork[J].Packaging and Food Machinery, 2007,25(5):53-55.

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XU Shao-hu,CUI Shuang.Analysis of Influencing Factors on Fruits MAP Inner Environment[J].Packaging Engineering, 2008,29(4):68-69.

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WANG Yuan-ming,ZHAO Jiang.The Reason on Changing Gas Composition and Detected Method on Gas in Package[J].China Packaging,2009(2):54-55.

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FAN Jun.The Relation between Headspace Gas Analysis and Shelf Life of Package[J].Packaging Engineering, 2012,33(5):146.

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WU Shu-yi,LIANG Yi.Observation and Discussion of the Latest In-line Headspace Oxygen Monitoring[J].Mechanical and Electrical Information,2011,(32):42-44.

[8] Jennifer M AMES,S Craig DUCKHAM,Jokie BAKKER.Choice and Use of Standards for Dynamic Headspace Trapping and Application to the Analysis of the Volatiles of Baked Potato[J].Headspace Analysis of Foods and Flavors, 2001,488:133-141.